


Tracing the fate of mercury from
vaccines in organisms:
methods for the determination of
ethylmercury, methylmercury, and
total mercury in biological tissues

Presenting Author: Annie Carter

Brooks Rand Labs

Seattle, WA, USA



Goal of the Study

Develop a procedure for the determination of MeHg, EtHg and Inorganic Hg in a tissue sample that yields acceptable recovery of SRMs and matrix spikes and yields reproducible results.

Thimerosal

- Vaccine preservative
- 50% (*m/m*) EtHg
- Been in use since 1930s
- In 1999, FDA recommended discontinuing use
- 2013 AAP recommends continuing the use

MeHg Analysis

- EPA 1630
- Used Brooks Rand Instruments MERX-M autoanalyzer



Sodium Tetraethylborate


“Ethylation” with NaBEt_4

- $\text{Hg}^0 \rightarrow \text{Hg}^0$ (no change-already volatile)
- $\text{CH}_3\text{-Hg}^+ \rightarrow \text{CH}_3\text{-Hg-CH}_2\text{-CH}_3$ (methylethyl mercury)
- $\text{Hg}^{2+} \rightarrow \text{CH}_3\text{-CH}_2\text{-Hg-CH}_2\text{-CH}_3$ (diethyl mercury)

Why NaBEt_4 can't be used to determine EtHg

- $\text{CH}_3\text{-CH}_2\text{-Hg}^+ \rightarrow \text{CH}_3\text{-CH}_2\text{-Hg-CH}_2\text{-CH}_3$ (diethyl mercury)

Chromatogram



Sodium Tetra(n-propyl)borate


“Propylation” with NaBPr₄

- Hg⁰ → Hg⁰ (no change-already volatile)
- CH₃-Hg⁺ → CH₃-Hg -CH₂-CH₂-CH₃ (methylpropyl mercury)
- Hg²⁺ → CH₃-CH₂-CH₂- Hg -CH₂-CH₂-CH₃ (dipropyl mercury)

Why NaBPr₄ can be used to determine EtHg

- CH₃-CH₂- Hg⁺ → CH₃-CH₂- Hg -CH₂-CH₂-CH₃ (ethylpropyl mercury)


Chromatogram




Propylation Reagent

- Purchased from ABCR
- Solution of 1g NaBPr₄, 2 g KOH, 97 mL DIW
- Aliquoted into 40mL amber glass vials (~10 mL/vial)
- Stored frozen at -20°C

Volume to Use



Propylation Reagent Holding Time



Correction Factor

- EPA Method 1630 section 12.4.2


$$F=100/R$$

F=Empirically derived correction factor


R=Recovery (the running mean of the last 30 IPR and OPR samples)

- Propylation is less efficient than ethylation
 - 1.15 (ethylated) and 1.39 (propylated) for MeHg
 - 1.49 for EtHg

Propylation Reagent Holding Time Study



Species Conversion



Sample Preparation – Experiment 1 TMAH

TMAH Digestion

- Add ~0.3 g sample and 2 mL 25% tetramethylammonium hydroxide (TMAH) solution to a glass vial. Bring volume to 10 mL with DI water.
- Seal the vial and heat to 65 °C in an oven for 4 hours. Cool to room temperature.
- Bring volume of samples up to constant level with DI water, to correct for volume loss.

Sample Preparation – Experiment 2 DCM

DCM Extraction

- Add sample and 5% H_2SO_4 + 18% KBr solution to a digestion vessel with a 1:4 (wt:v) ratio. Add 1 M CuSO_4 solution in a 1:5 (v:v) ratio with the H_2SO_4 /KBr solution.
- Add dichloromethane in a 2:1 (wt:v) ratio with the H_2SO_4 /KBr solution.
- Shake for 1 hour and let stand overnight.
- Centrifuge at 3,000 RPM for 15 mins. Separate the organic phase into a Teflon[®] bottle and add DI water in a ~5:2 (v:v) ratio. Heat to 70 °C until dichloromethane evaporates.

Sample Preparation – Experiment 3

KOH/MeOH plus DCM

- KOH/MeOH Digestion
- Add ~0.1 g sample and 1 mL 25% KOH/Methanol solution (wt/v) to a Teflon[®] digestion vial.
- Seal the vial and heat to 65 °C in an oven for 4 hours. Remove from the oven and bring the volume to 2.5 mL with methanol.
- Do DCM extraction procedure as previously described.


Sample Preparation – Experiment 4

KOH/MeOH plus Distillation


- Follow the KOH/MeOH Digestion procedure
- Transfer the digestate to a distillation vial
- Distill following EPA Method 1630 for water samples.

Blank Spikes


Blank Spikes




SRM Recovery



Matrix Spikes



Detection Limit




Additional Considerations

- incomplete digestion of liver and kidney matrix in the DCM prep
 - Because there is no certified reference material for ethylmercury in a solid tissue matrix, we could not be certain that this method was extracting all ethylmercury from the solid samples.


Kidney Sample Analysis

- MeHg and EtHg (simultaneous)
 - KOH/MeOH digest and distill
 - Propylation
- MeHg only (for confirmation)
 - KOH/MeOH digest
 - ethylation
- Total Hg by 1631e Appendix
 - 70:30 Nitric:Sulfuric digest
 - CVAFS on BRI MERX-T
- Inorganic mercury
 - Determined by difference

Sample results



% of THg that is EtHg



Batch QC

	Ethylmercury	Methylmercury
Mean Matrix Spike Recovery	106.9%	81.2%
Mean Matrix Spike Duplicate Precision	6.0%	10.6%
Mean Sample Duplicate Precision	12.7%	Non-detect
Mean Recovery NIST 955c SRM	86.8%	84.0%

Conclusions

- Substitution of NaBPr_4 for NaBEt_4 is an effective way to extend analytical equipment designed for MeHg analysis to the analysis of EtHg as well.
- The digest/distill method is effective for the simultaneous determination of EtHg and MeHg in liver and kidney samples
 - Methods developed for the analysis of blood and other bodily fluids are not necessarily applicable to tissue analysis due to their inability to disperse solid particles.
- There is a need for an animal tissue reference material certified for ethylmercury

Thanks!

Joel Creswell



Becky Thorsness

