



MINI

“Minimize Imperfections Naturally Improve”

\$299.00



Skin Care System utilizing sub-sensory current in micro amps. Multiple frequencies

Let your skin care professional establish a home skin care regimen incorporating the Mini. The Mini can be used for a deep cleansing, penetrate products for enhancement, and keep the skin in good condition, in between your professional treatments. Your skin care professional can suggest the correct products for your skin type and assist you in using the Mini for optimum results.

The Mini has an output of 330 μ A and uses Clareblend's multiple random frequency modulation for improved tone. The Mini has a constant current source to maintain the current output as well as safety features. The Mini runs on low dc voltage circuit for delivery of smooth current.

WARNINGS and Precautions:

The MINI facial procedure is intended for facial toning and should not be used if the following indications are present:

- Pregnancy
- Seizures
- Cancer or tumors
- Cardiac pacemaker, any electronic implanted defibrillators or stimulators or other electronic devices. (Possible interference with these devices may occur)
- Do not use if prone to phlebitis, thrombosis, or on varicose veins.

Use MINI only on normal skin. Do not use on skin that is:

- Swollen
- Infected
- Cancerous lesions.

This small hand held device, based on our professional system, provides quality and effectiveness.

Low current benefits have shown to:

- Increase circulation for skin health
- Increase Collagen and Elastin
- Minimize the appearance of fine lines and wrinkles
- Lift and firm the skin
- Enhance products with deeper penetration

Easy to use – “Plug and Go”

- No battery to replace
- Plugs into USB ports such as External chargers, computers, airport charging systems, or AC outlets (adapter included).
- Proven pre-determined settings provide great results without the guesswork!
- Gold connector for reliable power
- 1 year limited warranty

