




# MINI Testimonials



## Jodee B. NV 9/2014

I am very impressed with the Mini!! The size, the beep when it makes contact with the skin and then when the hold time is completed. The best thing though, is the results you achieve with this device. I expected great things from Clareblend but this has exceeded my expectations!!!

## Michon Kessler- NV - 10/1/14



Remove

Michon Kessler I love this thing!!!

October 1 at 7:52pm · Like

## Amy G. FL - 11/14/14

Love my mini- I actually use it at my desk, plugged into my computer as I read my morning email. I don't even have to think ,as the beep tells me what to do. Big timesaver as I can do my Microcurrent treatment at my desk. Great for a busy working woman like me.

## Michelle C. CA - 11/18/14

Remarkable Results Instantly

1. Tightens skin
2. Lifts eyes
3. Wrinkles around the eyes diminished
4. Looks like a Botox treatment
5. Makes skin feel good and alive
6. Face feels stimulated and radiant

I look forward to my treatments

**Carole W      CA - 11/11/14**

So convenient to use. The portability of the Mini makes it a tool that will be used often. I use the Mini for pain management for my husband's shoulder strains. He experiences results after just using it for 30 min. I use for facelift treatments. Love this product!

Sent from my iPhone

**Mary N.      OR 8/2014**

Very nice unit, easy to manage, easy instructions, good quality for price.

**Patricia T.      AZ      11/20/14**

As an Aesthetician, I use the Mini in most treatments. I always see immediate results, from product delivery to facial lifts. My clients love the results too! For personal use I love the benefits and it's easy to use.

**Linda W      CA      1/1/15**

I've had a full face lift (2 years ago) and wanted to keep the look going and my skin firm but fillers get costly. The Mini is perfect and it's my own natural stimulated collagen, Love it! P.S. My husband loves it on me too ☺ At 58, you can never look too young. LOL

**Cindi Benzer      CA      12/30/14**

I find it convenient and look forward to seeing the added improvement in addition to my regular appointments with Susan.

**Barb A.      OH      2/5/15**

Tried it last night and it was wonderful! It felt good and I actually noticed a difference in how my skin seemed to look.

**Jenill H.      GA      2/17/15**

Amazing results after 1 use. Thank You!

**Deborah Y.      AL      2/12/15**

Have been following Sunshine Botanicals, and was delighted to see the review of the Mini and said **YES!**

**Barbara C.      NC      2/23/15**

Mini seems to be reducing fine lines & wrinkles since I received it.

**Mary C.      CA      10-21-14**

I used the Clareblend Mini twice on my face & neck. A few days later at the salon my hairdresser said I looked younger and what did I do? The only thing different was the Mini treatment. You bet I'll keep using it. My husband used it on his sore knee, with a torn meniscus, he was in a lot of pain and the treatment helped reduce his pain. – Thank you

**Kim S.      CO      03/20/15**

Do you have a video on the mini? (Love it, BTW...)

**Lindsay      CA      3-29-16**

I have used it once since I received it yesterday. It was very easy to use and I already feel and see a difference. I love that it did not cause any discomfort. Excited to see more results !!

**Lovvet      CA      4-6-16**

I'm really enjoying it so far. I feel like my face looks fresh and lifted after each use. I am excited to see more results.

**Shirley Y      ON      3/18/16**

Happy that it visibly reduces the puffiness in the face. Helpful staff in fine tuning techniques.

**Linda D.      TN      6/15/16**

This product is truly amazing! I was skeptical, but I saw a difference in Christine's face in her video, and if you could see a change on such a pretty face, I thought I would try it. I don't do anything special for skin care so this was all new to me. I am shocked at how different my face looked just after 30 minutes of using it. It truly changed the structure of my face, and it even feels different – toned and tighter. I have my defined jaw line back. I am SO incredibly happy and blessed that I found this device!! Thank you for making such a great product!!

**Carmalee E.      UT      6/14/16**

I've only used it once on myself, but am very excited about the prospect of its results and what it will do!