

**HEARTSINE SAMARITAN PUBLIC ACCESS DEFIBRILLATOR 300/300P:**

**CLASS I RECALL**

HeartSine notified customers that certain Samaritan 300/300P PAD devices have been found to intermittently turn on and off, which may eventually deplete the battery. In addition and separately, certain Samaritan 300/300P PAD devices containing early versions of the battery management software may misinterpret a temporary drop in battery voltage as signaling a low battery, resulting in the device turning itself off. A device experiencing either condition could be unable to deliver therapy during a cardiac event. For more information on this recall go to the FDA web page at <http://www.fda.gov/Safety/MedWatch/SafetyInformation/SafetyAlertsforHumanMedicalProducts/ucm328614.htm>